


# DENVER NASAL SPLINT



*Stabilise nasal bones and cartilage  
and reduce oedema after nasal surgery*

The Denver range consists of **three** different types of splint to suit the demanding requirements for stabilisation following Rhinoplasty Surgery. By using a **Denver** splint you will save valuable time as application takes 60 seconds only. The original Denver splint is not to be confused by inferior imitations.

This 7-10 day splint combines adhesive layers with ventilating Velcro giving gentle pressure to the subcutaneous tissue, reducing oedema.

Easy to apply, the Denver Nasal Splint® is:

- **Inconspicuous**
- **Water resistant**
- **Fast, usually one minute application**
- **Clean, no messy application**
- **Easy to use kit**

For more information on the Denver Nasal Splint® range of products, please contact:

Phillip Silver, Associate Product Manager

---

**FORTH MEDICAL LIMITED**

Forth House • 42 Kingfisher Court • Hambridge Road Newbury • Berkshire • RG14 5SJ

Telephone: 01635 550100 • Fax: 01635 550050

Email: [phil@forthmedical.ltd.uk](mailto:phil@forthmedical.ltd.uk)

[www.forthmedical.com](http://www.forthmedical.com)

# DENVER NASAL SPLINT


## INDICATIONS:

All Denver Nasal Splints® are used in the stabilisation of the nose and the reduction of oedema after nasal surgery or trauma.

## BIOMECHANICS:

The moldable yet firm aluminum retains the nasal bones in position and the Velcro™ or Velfoam™ layers give cushioned pressure on subcutaneous tissue, reducing the Oedema. Packaged non-sterile with sterile components.

## KIT COMPONENTS



- Instruction sheet
- 1 x alcohol wipe
- 1 x skin prep
- 1 x dorsal pad
- Pre-cut Protecto™ Nose Strips
- 1 x Denver Adhesive Nasal Splint®

## DESIGNER COLOURS!



All Denver Splint® models are available in Designer colours. At no extra charge, designer splints are fun for teenagers, athletes, and rhinoplasty patients!

For more information on the Denver Nasal Splint range of products, please contact:

Phillip Silver, Associate Product Manager

## FORTH MEDICAL LIMITED

Forth House • 42 Kingfisher Court • Hambridge Road Newbury • Berkshire • RG14 5SJ

Telephone: 01635 550100 • Fax: 01635 550050

Email: [phil@forthmedical.ltd.uk](mailto:phil@forthmedical.ltd.uk)

[www.forthmedical.com](http://www.forthmedical.com)